

## EYEPIECE GLOBE CAM



### Globecam V (video)

Product reference. 9855.1019

The eyepiece camera can easily be inserted into the ocular tube of any microscope.

For convenience a parfocal image is displayed on screen.

### Digital Camera on Flexible Arm

Product Reference. 9544.7200




## Si-3000



3 MP low-noise digital camera with integrated capture & analysis software.

Product reference. 9556.3000

- High performance 3 Mpx camera development for microscopy
- 1/2" 2048 x 1536 (QXGA) image array and on-chip converter
- Up to 27 frames per second (fps)
- Full resolution and 40fps at XGA (1024 x 768) resolution
- Sensitivity 1.0V/Lux-sec
- Complete with USB cable and image analysis software
- Tripod fitting



## DIGITAL USB INSPECTION CAMERA



- High performance 3 Mpx camera mounted on sturdy adjustable stand
- 1/2" 2048 x 1536 (QXGA) image array and on-chip converter
- Up to 27 frames per second (fps)
- Full resolution and 40fps at XGA (1024 x 768) resolution
- Complete with USB cable and image analysis software
- LED ring illuminator high clarity 4x Objective
- Plug and play software
- Option of 4x objective or Pentax adjustable optic system