

DIGITAL SLR CAMERA



8MP Digital SLR camera.
Product reference. 9780.3000

Recommended accessory:
Product reference. 9780.2222
 1.6x EOS M42 'C' mount adapter with T2 ring

Model and specification subject to change without prior notification.

* Please remember to order the necessary C-mount adapter



IMAGE ANALYSIS PACKAGE



The ultimate software package for photo-microscopy.

Product reference. 9600.0011

- Capture & enhancement of images
- Full analysis & measuring tools/techniques

Image Pro Express is a cost effective software package perfect for routine imaging. The software comes equipped with all the features needed to capture and enhance images for scientific, medical and industrial research.

Image processing, measurement and analysis software.

Exceptional levels of software technology and analysis tools will allow you to streamline image capture, measurements and enhancements whilst improving the accuracy of your results.

All packages feature (as standard):

Live Measurement and Overlay Settings / Calibration (automatic & manual) / Extended Focus Imaging (EFI) / 3D Visualisation / Image Stitching / Live Image Comparisons / Time Lapse Capture and Movie File Production / Combine Image Planes (Fluorescence Imaging) / Import data into Excel / Large selection of Measuring and Analysis Functions (including annotation).

Available in 3 Packages:

i-Solution Lite: Manual measurements, including moderate image processing functions.

Product reference. 9600.1010

i-Solution: Advanced image processing and enhancement, including both manual and automated image analysis.

Product reference. 9600.1012

i-Solution DT: Premium software package incorporating all elements of image analysis found in i-Solution and i-Solution Lite, including motorized stage control and 3D analysis. Customized modules can also be added by request.

Product reference. 9600.1013

